

Safety of orthopedic implants: implant migration analysis a must

Hasan, S.

Citation

Hasan, S. (2024, June 4). *Safety of orthopedic implants: implant migration analysis a must*. Retrieved from https://hdl.handle.net/1887/3762018

Version:	Publisher's Version			
License:	<u>Licence agreement concerning inclusion of doctoral</u> <u>thesis in the Institutional Repository of the University</u> <u>of Leiden</u>			
Downloaded from:	https://hdl.handle.net/1887/3762018			

Note: To cite this publication please use the final published version (if applicable).



Chapter II

TKA implants on average have lower mean 10-year revision rates than non-RSA-tested designs

> Shaho Hasan¹, Perla J Marang-van de Mheen², Bart L Kaptein¹, Rob GHH Nelissen¹, Bart G Pijls¹

¹Department of Orthopaedic, Leiden University Medical Center, Leiden, The Netherlands ²Department of Biomedical Data Sciences, Leiden University Medical Center, Leiden, The Netherlands

> Clin Orthop Relat Res. 2020 Jun;478(6):1232-1241. doi: 10.1097/CORR.000000000001209

Abstract

Background

The number of revisions after TKA is expected to rise because of aging populations in many countries and because patients are undergoing TKA at younger ages. Aseptic loosening is a major reason for late revision, which can be predicted by radiostereometric analysis (RSA) of small groups of patients at 2 years of follow-up. RSA is therefore an ideal tool to assess new TKA designs before they are introduced to the market, although not every TKA design has been studied with RSA. If RSAtested TKA designs have lower 10-year revision rates in national registries than non-RSA-tested TKA designs, RSA testing of all new designs could be advocated.

Questions/purposes

In this study, we asked: Is there a difference in the all-cause revision rate between non-RSA-tested and RSA-tested TKA designs registered in national knee arthroplasty registries at 5 and 10 years of follow-up?

Methods

Knee arthroplasty registries were identified through the European Federation of National Associations of Orthopaedics and Traumatology webpage and through a manual internet search. Inclusion criteria were a minimum follow-up duration of 10 years and available revision or survival data per TKA design. Twenty-six registries were identified; seven were included comprising 339 TKA designs, of which 236 designs were classified as RSA-tested and 103 as non-RSA-tested. Six registries were excluded because no report was published. One registry was excluded because no fixation method was mentioned (79 TKA designs). Another registry was excluded because there was no 10-year data available (22 non-RSA-tested designs; 10 RSAtested designs). Eleven registries were excluded because these registries did not provide revision rates per design and had not reached 10 years follow-up. The revision rates with their standard errors were extracted per design. We used the data from a recent meta-analysis to identify whether a TKA design was previously tested with RSA. This meta-analysis found 53 RSA studies comprising 70 different TKA designs. The prosthesis model, fixation method and insert type were extracted from these RSA-studies. The design characteristics of the TKA reported in the knee arthroplasty registries were also extracted, and if possible, matched to the TKA designs reported in the RSA-studies. At 5 years of follow-up, 191 TKA designs were identified as non-RSA-tested and 92 were identified as RSA-tested. At 10 years of follow-up, 154 TKA designs and 74 TKA designs were classified as non-RSA-tested and RSA-tested, respectively. A random-effects model using the Metafor Package in R statistics was used to estimate the pooled revision rate at 5 and 10 years of followup for both groups. The difference in revision rates between groups at 5 and 10 years of follow-up was estimated by including RSA as a factor in the random-effects model.

Results

Mean all-cause revision rates at 5 years for non-RSA-tested and RSA-tested implants were 3.6% (95% CI 3.4 to 3.8) and 2.9% (95% CI 2.7 to 3.0), with a mean difference of 0.6% favoring RSA-tested implants (95% CI 0.4 to 0.8; p < 0.001). Mean all-cause revision rates at 10 years for non-RSA-tested and RSA-tested implants were 5.5% (95% CI 5.2 to 5.9) and 4.4% (95% CI 4.1 to 4.7), with a mean difference of 0.9% favoring RSA-tested implants (95% CI 0.4 to 1.3; p < 0.001).

Conclusions

Although there are exceptions, across registries, TKA designs that have been tested in an RSA setting have a slightly lower (about 1%) mean all-cause revision rate at 5year and 10-year follow-up than those tested in a non-RSA setting do. Acknowledging the inherent limitations of this observational study, a risk difference of 1% could potentially translate into an approximate 20% decrease in revision burden up to 10 years, which may have a profound impact on patient morbidity and health-related costs.

Level of Evidence

Level III, therapeutic study.

Introduction

The number of revisions after TKA is expected to rise because of aging populations in many countries, and because of increased usage of this procedure in younger patients.^{1, 2} Unfortunately, the introduction of newer TKA designs has not always resulted in fewer revisions.^{3, 4} A major reason for long-term revision of a TKA implant is aseptic loosening, which can be predicted using the 2-year postoperative prosthesis migration profile, measured using radiostereometric analysis (RSA).^{5, 6} RSA was first described in 1974, has been improved for use with digital radiographs, and has been used with various TKA designs.⁷⁻¹³ Given the high precision of RSA, RSA studies generally need only approximately 50 patients per group to detect a difference in prosthesis migration between TKA designs, making RSA an ideal tool to evaluate new TKA designs in early clinical trials.¹⁴ The importance of RSA studies before widespread market introduction of new designs has been noted in numerous reports that correlate early (1 to 2 years) migration patterns of knee implants with 10year survival of these implants.^{5, 6, 14, 15} Phased introduction of new TKA designs, including those evaluated in early clinical RSA trials, has been proposed to improve patient safety.14, 16-19

However, not every TKA design has been studied with RSA before market introduction. In the AOANJR registry for instance, nearly 194 different TKA design combinations have been registered, with reported 10-year survival rates ranging from 86.5% to 98.1%, and most designs were not evaluated in an RSA study.²⁰ RSA could be used to warn clinicians about implants that are more likely to have an increased risk of aseptic loosening, thus safeguarding against the widespread use of such implants. Such a warning might result in withdrawal of designs from the market, thereby leaving only the better-performing implants and preventing many early revisions.¹⁴ Following this mechanism, TKA designs tested with RSA may be expected to have a lower revision rate during long-term follow-up than non-RSAtested TKA designs. In an earlier report with shorter follow-up, Nelissen et al. (2011) found that RSA-tested TKA designs had a lower revision rate in three national knee arthroplasty registries with up to 5 years of follow-up.14

Here, we used six national registries and one regional registry to answer the question: Is there a difference in the all-cause revision rate between non-RSA-tested and RSA-tested TKA designs registered in national knee arthroplasty registries at 5 and 10 years of follow-up?

Materials and Methods

Study Search

Through the Network Orthopaedic Registries of Europe—European Federation of National Associations of Orthopaedics and Traumatology webpage (EFORT), national and regional knee arthroplasty registries were identified.²¹ A manual worldwide-web search was then conducted to identify any knee arthroplasty registry not listed on the EFORT webpage. Published reports were extracted from these registries. Inclusion criteria were a minimum follow-up duration of 10 years and available revision or survival data for each TKA design. Knee arthroplasty registries were excluded if no information regarding the fixation method was provided. However, if a study or report stated that more than 90% of the TKA designs were cemented, the entire registry was included and all TKA designs were assumed to be cemented (but tested in a sensitivity analysis – see below). No language restriction was used.

The search yielded 26 annual reports of knee arthroplasty registries, of which six national registries (from Australia, Finland, New Zealand, Norway, Sw eden, and the United Kingdom) and one regional registry (Emilia-Romagna, Italy) were included (Fig. II.1).^{20, 22-27} TKA designs of one registry were excluded due to unknown fixation method. TKA designs from another registry were excluded as only 7-year data was available. All other excluded registries did not have 10-year follow-up and did not report revision rates per TKA design (Fig. II.1). From the seven registries, 339 TKA designs were extracted. The maximum follow-up duration ranged between 13 and 41 years, and all registries had a completeness of \geq 95% for primary TKA. The definition

of completeness was not clarified in all registries but was defined as the percentage of patients receiving a primary TKA included in the registry in most registries. The mean age at the time of surgery ranged from 68 years to 71 years [Table II.I]. The Finnish registry did not report a mean age but divided patients into four age groups (younger than 55 years, 55 to 64 years, 65 to 74 years, and 75 years and older), with most patients (39%) in the 65 to 74 years age group. The proportion of female patients in the registries ranged from 47% in the Sweden registry to 71% in the Emilia-Romagna, Italy registry [Table II.I].

		Australia	Finland	Emila-	New	Norway	Sweden	United
				Romagna	Zealand			Kingdom
				(Italy)				
TKAs (n)		547,407	194,787	39,782	93,497	29,834	109,393	975,739
Follow-up (year	rs)	16	25	16	17	23	41	13
TKA designs (n	l)	143	356	39	34	19	13	56
Publication yea	r	2017	2018	2017	2017	2018	2017	2017
Completeness (%)	98%	96%	98%	>95% ^a	97%	97%	96%
Age (mean, yea	rs)	68.5	65-74	70.6	68	68.5	69	70
Sex (female, %))	56%	68%	71%	52%	63%	47%	57%
The three most-used TKA								
designs								
	1	Triathlon ^b	Triathlon ^b	Attune ^d	Triathlon ^b	NexGen ^c	NexGen ^c	
	2	NexGen	NexGen ^c	NexGen ^c	Attune ^d	LCS	PFC ^d	
		Flex CR ^c				Completed		
	3	NexGen	PFC	Legione	Genesis II ^e	PFC	Triathlon ^b	
		Flex LPS ^c	Sigma ^d			Sigma ^d		
Revision due to	loosening of	26%	9%		19%		26%	26%
all TKA in regi	stry (%)							
TKAs (n)		547,407	194,787	39,782	93,497	29,834	109,393	975,739

Table II.I Characteristics of included registries

aIn 95% of public hospitals.

^bStryker Inc, Mahwah, NJ, USA

Cimmer Inc, Warsaw, IN, USA

^dDepuy Synthes, Warsaw, IN, USA

°Smith&Nephew, Memphis, TN, USA

reconstruction prostheses were excluded. The number of designs ranged between 13 and 143 per annual report.

To identify whether a TKA design was previously tested with RSA, we used the data from a recent meta-analysis.¹⁵ In short, this meta-analysis searched PubMed, EMBASE, Web of Science, and the Cochrane Library for studies using RSA and primary TKA before July 2016. Data on all designs were extracted from the 53 included studies, which included 70 different RSA-tested TKA designs [Fig. II.1]. Design characteristics of the TKA reported in the knee arthroplasty registries were extracted and, if possible, matched to the TKA designs reported in the RSA-studies. Every TKA design reported in the included knee arthroplasty registries was classified as non-RSA-tested or RSA-tested, resulting in two groups in every registry. For a design to be classified as RSA-tested, the design in the registry had to be identical to the design reported in an RSA study. If the insert was not specified in the registry, but the design and fixation matched the TKA design, the design was classified as RSA-tested [Fig. II.1].

Seven registries with 339 TKA designs were included of which 236 were classified as non-RSA-tested and 103 as RSA-tested. Fixation was uncemented for 54 designs and cemented for 285 designs. Cruciate-retaining inserts were used in 110 designs and posterior-stabilizing in 72 designs or not explicitly mentioned. Mobile bearings were used in 49 designs. The Norwegian registry only reported 3-year and 10-year revision rates and could therefore not be included at 5 years. At 5 years, 191 TKA designs were identified as non-RSA-tested [Supplemental data II.I]. In addition, 92 TKA designs were identified as RSA-tested [Supplemental data II.I]. At 10 years, 154 TKA designs were identified as non-RSA-tested [Supplemental data II.II] and 74 designs were identified as RSA-tested [Supplemental data II.IV]. Between baseline and 10 years, 82 non-RSA-tested and 29 RSA-tested designs could not be included as these had not reached 10 years of follow-up.

Figure II.I Inclusion Flowchart



Data Analysis

First, a random-effects model was used to calculate the pooled all-cause revision percentages and their standard errors at 5 and 10 years of follow-up for the non-RSA-tested and RSA-tested TKA designs, including a DerSimonian-Lard estimator to take into account the heterogeneity between the designs.²⁹ RSA-tested (yes or no) was included as a factor to test for a difference between groups at 5 and 10 years of follow-up. Moreover, pooled all-cause revision percentages and their standard errors were calculated separately for each registry for non-RSA-tested and RSA-tested TKA designs.

In all random-effects models, a DerSimonian-Lard estimator was used to estimate heterogeneity.²⁹ In a sensitivity analysis, the more conservative empirical Bayes estimator was used to test whether the heterogeneity estimator would affect the results.³⁰ The I² was used to estimate the extent of heterogeneity. The heterogeneity is the variation between the designs in both groups, which is considered low, moderate, or high if I² is 25%, 50% or 75%, respectively.^{31, 32} Outcomes are given in percentages with 95% CIs. The Metafor Package in R Statistics (version 3.6.1; R Foundation for Statistical Computing, Vienna, Austria) was used for all analyses.³³

Post-hoc sensitivity analyses

We performed three post-hoc sensitivity analyses to test the impact of various assumptions on the primary outcome. The first analysis excluded data from registries for which the fixation method was missing (Sweden and Emilia-Romagna, Italy), and data from the registry that did not report the insert of the design (New Zealand). The second analysis included the four RSA studies from the meta-analysis that were excluded from the primary analysis because of not reporting migration data or other reasons. This resulted in reclassification of six non-RSA-tested TKA designs as RSA-tested. The third analysis included data from the Danish and Dutch knee arthroplasty registries that fulfilled all but one of the inclusion criteria. ^{34, 35} The Danish TKA designs lacked information on fixation and were assumed to be

cemented in this sensitivity analysis, and the Dutch registry published 10 years follow-up data in November 2019 (after initial manuscript submission) and could only be included recently.

Results

Revision Rates of Non-RSA-tested and RSA-tested TKA Designs at 5 Years of Follow-up

All-cause revision at 5 years was slightly less for the RSA-tested designs than for the non-RSA-tested designs [Fig. II.II]. Mean all-cause revision rates at 5 years for non-RSA-tested and RSA-tested implants were 3.6% (95% CI 3.4 to 3.8) and 2.9% (95% CI 2.7 to 3.0), with a mean difference of 0.6% (95% CI 0.4 to 0.8; p < 0.001) favoring RSA-tested implants. Using the more conservative Empirical Bayes estimator, the mean difference was 0.7% (95% CI 0.3 to 1.0; p < 0.001) in favor of RSA-tested implants.

The revision rates of the RSA-tested TKA designs in the registries ranged between 2.3% and 3.9%, whereas revision rates of the non-RSA-tested TKA designs ranged from 2.5% to 4.7%. In all registries, the point estimate of RSA-tested TKA designs was lower than that of non-RSA-tested designs, but the absolute difference between groups was smallest in the United Kingdom (0.2% at 5 years of follow-up). The highest revision rate of RSA-tested implants was reported in Finland (3.9% at 5 years of follow-up). New Zealand and Sweden used more RSA-tested TKA designs than non-RSA-tested TKA designs, in contrast to other countries. Australia had the greatest number of TKA designs registered (n = 126). Within the RSA-tested and non-RSA-tested groups, high variation was found between the TKA designs, expressed by the high heterogeneity (I² 96% in the non-RSA-tested group and 98% in the RSA-tested group). Including fixation or insert in the model did not reduce the heterogeneity, suggesting that there is large variation in revision rates between designs. In addition, it is important to note that although a slightly lower mean all-cause revision rate was found for RSA-tested TKA, some non-RSA-tested TKA

performed well, whereas some RSA-tested TKA performed poorly.

Fig. II.II This forest plot shows revision rates of the non-RSA-tested and RSA-tested TKA designs with 95% CIs at 5 years of follow-up per registry and the pooled revision rate per group with 95% CI. In addition, the sensitivity analysis including designs from the Dutch and Danish knee arthroplasty registry.



Revision Rates of Non-RSA-tested and RSA-tested TKA Designs at 10 Years of Follow-up

Similarly, all-cause revision at 10 years was slightly less common among RSA-tested designs than it was in non-RSA-tested designs [Fig II.III]. Mean all-cause revision rates at 10 years for non-RSA-tested and RSA-tested implants were 5.5% (95% CI 5.2 to 5.9) and 4.4% (95% CI 4.1 to 4.7), with a mean difference of 0.9% (95% CI 0.4 to

1.3; p < 0.001) favoring RSA-tested implants. Using the more conservative Empirical Bayes estimator, the mean difference was 0.9% (95% CI 0.2 to 1.6; p = 0.01) favoring RSA-tested implants. The revision rates in the registries ranged between 3.9% and 8.0% for non-RSA-tested and between 3.6% and 6.4% for RSA-tested TKA designs with large heterogeneity in both groups (I^2 97%).

Fig. II.III This forest plot shows revision rates of the non-RSA-tested and RSA-tested TKA designs with 95% CIs at 10 years of follow-up per registry and the pooled revision rate per group with 95% CI. In addition, the sensitivity analysis including designs from the Dutch and Danish knee arthroplasty registry.



Revision percentage (%) with 95% Confidence Interval

Post-hoc Sensitivity Analyses

First, excluding the data from registries with assumed fixation method or inserts (Sweden, New Zealand and Emilia-Romagna, Italy) resulted in a slightly smaller mean difference in all-cause revision rate between groups of 0.5% (95% CI 0.2 to 0.8; p < 0.001) at 5 years and 0.7% (95% CI 0.2 to 1.2; p = 0.003) at 10 years. Second, reclassifying the TKA-designs from the excluded studies did not influence the mean revision rates in both groups nor on the difference between groups (data not show n). Third, including both the Danish and Dutch registries, the mean difference of all-cause revision rate between RSA-tested and non-RSA-tested designs was 0.6% (95% CI 0.4 to 0.8; p < 0.001) in favor of RSA-tested designs at 5-year follow-up. At 10-year follow-up, the mean difference in all-cause revision was 0.9% (CI 95% 0.4 to 1.3; p < 0.001) favoring RSA-tested implants.

Discussion

Regulations regarding the introduction of new orthopaedic devices should have a healthy balance between innovation and patient safety.¹⁸ To improve patient safety, new medical device regulations were established in Europe; they require clinical evidence before new implants are introduced to the European Union market.³⁶ RSA may be an important part of such clinical testing, and its use as an early-warning system for implants likely to fail as a result of aseptic loosening has often been proposed.^{14, 16, 18, 37, 38} However, it is unknown whether RSA-tested TKA designs are associated with a lower revision rate during long-term follow-up in registries, though this may seem likely if problematic RSA tested designs are withdrawn from the market. By pooling data from several national registries and a regional registry, we found that implants that had undergone RSA testing, overall had a slightly (about 1%) lower all-cause revision rate at 5 and 10 years compared with implants that had not undergone RSA testing.

We should consider the following limitations. First, our study was an observational study and cannot imply causation between RSA and a lower TKA revision rate, but

rather showed an association between these two factors. Second, the classification of TKA designs as RSA-tested or non-RSA-tested came from another meta-analysis.¹⁵ However, a post-hoc sensitivity analysis showed similar results after reanalyzing the data from the meta-analysis and reclassifying the six TKA-designs that were excluded in the meta-analysis from non-RSA-tested to RSA-tested. Third, we should consider the possibility that differential loss to follow-up may have influenced the results, although here the loss to follow-up was comparable between both groups. Fourth, revision rates as reported from the registries were used as the outcome measure, which is relatively crude. Such rates are influenced not only by the performance of a particular TKA design, but also by patients' complaints (for example, pain) and the surgical decision-making process, which is affected by factors such as patients' comorbidities, cultural differences between patients (such as pain acceptance), and waiting lists.³⁹ Nevertheless, many implants were included in the study, and we assume that both groups were similarly affected by these factors influencing revision. Fifth, mechanical loosening of the tibia is not the only reason for revision. Other common reasons are instability and infection, which are not assessed by RSA.35 A phased introduction should therefore include clinical trials to assess these contributing factors for revision. In addition, the absolute difference was small at both 5 and 10 years (less than 1%), raising the question of the relevance of this effect. However, this effect should be interpreted considering the total revision rate, which is also low (about 5% at 10 years), meaning an absolute difference of 0.5% to 1% results in a decrease of approximately 10% to 20% for all-cause revision at 5 and 10 years. Considering the enormous number of TKA procedures performed globally, a 1% decrease in TKA revision could have a tremendous impact on the burden for patients needing a TKA revision and result in considerable reduction of health-related costs. Another limitation that should be noted is that there was high heterogeneity in all analyses, which could not be explained by the fixation method or the different inserts (data not shown). Heterogeneity is thus likely attributed to the many different designs included in the study with varying performance between the different designs. It should thus be emphasized that not all non-RSA-tested TKA

Chapter II

designs performed poorly and, vice versa, not all RSA-tested TKA designs performed well. Finally, we had to assume the fixation method for two registries and the insert type for one registry, which might have results in misclassification of some TKA designs although our sensitivity analysis showed this was not likely to change the results or conclusions.

We found that RSA-tested TKA designs had a slightly (about 1%) lower mean allcause revision rate at 5 and 10 years than non-RSA-tested designs. These results are in line with a previously published study comparing non-RSA-tested and RSA-tested TKA designs in three knee arthroplasty registries up to 5-year follow-up.¹⁴ Our findings might be explained by the fact that RSA could provide an early warning about inferior TKA designs that fail because of aseptic loosening of the tibia. This early warning function could theoretically lower revision rates if poorly performing implants were withdrawn from the market or no longer used, and well-studied and excellent-performing TKA designs continuing to be used. Given our observational study we were unable to test this hypothesis in the present study or determine whether this is the case. Possible alternative explanations could be that RSA testing is a proxy for a rigorous clinical testing program by the manufacturer, or that more prudent surgeons are more likely to use RSA tested implants.

Before introduction of the new European medical device regulations, a phased introduction of new implants was proposed by several authors and the Idea, Development, Exploration, Assessment, Long-term Study—Devices (IDEAL) consortium to guide the introduction of novel devices.^{14, 16, 18, 37, 38, 40} The best clinical introduction of a new TKA implant, in our opinion, would be to clinically evaluate implant fixation (that is, micromotion) as well as the surgical procedure. Thus, RSA studies and larger prospective studies could be nested in national or regional registries. Beyond Compliance, an initiative originating from the United Kingdom supporting the safe introduction of implants by bringing clinicians, implant manufacturers and an independent expert panel together to assess outcomes of joint replacements, could be performed parallel to RSA studies.^{18, 41} RSA studies or implant migration studies (Einzel-Bild-Röntgen-Analyse, CT-RSA) could play an important role in such a phased, stepwise introduction of new implants.⁴² Because of the accuracy of the RSA technique, there is no need to expose large groups of patients to new implant designs that could potentially be inferior to the current state-of-the-art designs. In addition to exposing fewer patients, shorter follow-up is needed, as migration results of implants after 2 years are often able to show differences in migration, in contrast to the long-term follow-up needed for classic observational studies, with survival of the implant as endpoint.^{14, 16}

Reducing the revision rate of TKA is particularly of interest because this procedure is estimated to increase by approximately 600% between 2005 and 2030, resulting in 268,200 revisions in 2030 in the United States alone.¹ Considering that the mean cost of revision TKA in the United States is USD 49,360, using only selected wellperforming TKA designs might save billions of dollars annually.⁴³

Conclusions

The number of different TKA designs is enormous, with new designs being introduced almost annually, and surgeons should remain skeptical about novel designs without proper evidence.^{4, 20} Several well-studied and excellent-performing TKA designs are currently available, and new designs should prove that they outperform these legacy products before replacing them. RSA testing is one method of testing new prosthesis introductions. Although there are exceptions, we found that TKA designs tested in an RSA setting were associated with a slightly lower (about 1%) mean all-cause revision rate at 5-year and 10-year follow-up than those tested in a non-RSA setting. The relevance of this small effect should be interpreted in the context of this being a relative decrease of approximately 20% for all-cause revision at 5 and 10 years while also considering the enormous number of TKA procedures performed globally. Future studies should address the possible explanations for the association found between RSA-testing and a lower mean all-cause revision.

References

- Kurtz S, Ong K, Lau E, Mowat F, Halpern M. Projections of primary and revision hip and knee arthroplasty in the United States from 2005 to 2030. J Bone Joint Surg Am. 2007 89(4):780-5 DOI: 10.2106/jbjs.f.00222.
- Patel A, Pavlou G, Mujica-Mota RE, Toms AD. The epidemiology of revision total knee and hip arthroplasty in England and Wales: a comparative analysis with projections for the United States. A study using the National Joint Registry dataset. Bone Joint J. 2015 97(8):1076-81 DOI: 10.1302/0301-620X.97b8.35170.
- Anand R, Graves SE, de Steiger RN, Davidson DC, Ryan P, Miller LN, et al. What is the benefit of introducing new hip and knee prostheses? J Bone Joint Surg Am. 2011 93:51-4 DOI: 10.2106/jbjs.K.00867.
- Nieuwenhuijse MJ, Nelissen RG, Schoones JW, Sedrakyan A. Appraisal of evidence base for introduction of new implants in hip and knee replacement: a systematic review of five widely used device technologies. Bmj. 2014 349:85133 DOI: 10.1136/bmj.g5133.
- Ryd L, Albrektsson BE, Carlsson L, Dansgard F, Herberts P, Lindstrand A, et al. Roentgen stereophotogrammetric analysis as a predictor of mechanical loosening of knee prostheses. J Bone Joint Surg Br. 1995 77(3):377-83.
- Pijls BG, Valstar ER, Nouta KA, Plevier JW, Fiocco M, Middeldorp S, et al. Early migration of tibial components is associated with late revision: a systematic review and meta-analysis of 21,000 knee arthroplasties. Acta Orthop. 2012 83(6):614-24 DOI: 10.3109/17453674.2012.747052.
- Selvik G. Roentgen stereophotogrammetry. A method for the study of the kinematics of the skeletal system. Acta Orthop Scand Suppl. 1989 232:1-51.
- Kaptein BL, Valstar ER, Stoel BC, Rozing PM, Reiber JH. A new model-based RSA method validated using CAD models and models from reversed engineering. J Biomech. 2003 36(6):873-82.
- 9. Valstar ER, Vrooman HA, Toksvig-Larsen S, Ryd L, Nelissen RG. Digital automated

RSA compared to manually operated RSA. J Biomech. 2000 33(12):1593-9.

 Molt M, Toksvig-Larsen S. 2-year followup report on micromotion of a short tibia stem. A prospective, randomized RSA study of 59 patients. Acta Orthop. 2015 86(5):594-8 DOI: 10.3109/17453674.2015.1033303.

11.

Carlsson A, Bjorkman A, Besjakov J, Onsten I. Cemented tibial component fixation performs better than cementless fixation: a randomized radiostereometric study comparing porous-coated, hydroxyapatite-coated and cemented tibial components over 5 years. Acta Orthop. 2005 76(3):362-9.

 Laende EK, Richardson CG, Dunbar MJ. A randomized controlled trial of tibial component migration with kinematic alignment using patient-specific instrumentation versus mechanical alignment using computer-assisted surgery in total knee arthroplasty. Bone Joint J. 2019 101(8):929-940 DOI: 10.1302/0301-620X.101b8.Bjj-2018-0755.R3.

 Tjornild M, Soballe K, Hansen PM, Holm C, Stilling M. Mobile- vs. fixed-bearing total knee replacement. Acta Orthop. 2015 86(2):208-14 DOI: 10.3109/17453674.2014.968476.

 Nelissen RG, Pijls BG, Karrholm J, Malchau H, Nieuwenhuijse MJ, Valstar ER. RSA and registries: the quest for phased introduction of new implants. J Bone Joint Surg Am. 201 93:62-5 DOI: 10.2106/ibis.k.00907.

 Pijls BG, Plevier JWM, Nelissen R. RSA migration of total knee replacements. Acta Orthop. 2018:1-9 DOI: 10.1080/17453674.2018.1443635.

- Malchau H. Introducing new technology: a stepwise algorithm. Spine (Phila Pa 1976). 2000 25(3):285 DOI: 10.1097/00007632-200002010-00004.
- 17. Malchau H, Garellick G, Berry D, Harris WH, Robertson O, Karrlholm J, et al. Arthroplasty implant registries over the past five decades: Development, current, and future impact. J Orthop Res. 2018 36(9):2319-2330 DOI: 10.1002/jor.24014.
 18. Sedrakyan A, Campbell B, Merino JG,

Kuntz R, Hirst A, McCulloch P. IDEAL-D:

	a rational framework for evaluating and		https://ripo.cineca.it/authzssl/pdf/report
	regulating the use of medical devices. Bmj.		eng 2016.pdf. Accessed on July 17, 2018.
	2016 353;i2372 DOI: 10.1136/bmi.i2372.	28.	Hozo SP, Diulbegovic B, Hozo I.
19.	Huiskes R. Failed innovation in total hip		Estimating the mean and variance from
	replacement. Diagnosis and proposals for		the median, range, and the size of a
	a cure. Acta Orthop Scand. 1993		sample. BMC Med Res Methodol. 2005 5:13
	64(6):699-716 DOI:		DOI: 10.1186/1471-2288-5-13.
	10.3109/17453679308994602.	29.	DerSimonian R. Laird N. Meta-analysis in
20.	Australian Orthopaedic Association.		clinical trials. Control Clin Trials. 1986
	National Joint Replacement Registry: Hip.		7(3):177-88 DOI: 10.1016/0107-
	Knee & Shoulder Arthroplasty Annual		2456(86)90046-2.
	Report 2017. Available from:	30.	Morris CN. Parametric Emprical Bayes
	https://aoanirr.sahmri.com/documents/101	,	Inference: Theory and Applications.
	80/397736/Hip%2C%20Knee%20%26%20S		Journal of American Statistical
	houlder%20Arthroplasty. Accessed on		Association. 1983 78:381:47-55.
	July 20, 2018.	31.	Higgins JP, Thompson SG. Quantifying
21.	Network of Orthopaedic Registries of		heterogeneity in a meta-analysis. Stat
	Europe, NORE research & publications.		Med. 2002 21(11):1539-58 DOI:
	Available from:		10.1002/sim.1186.
	https://www.efort.org/about-	32.	Higgins JP, Thompson SG, Deeks JJ,
	us/nore/research/. Accessed on July 14.	,	Altman DG. Measuring inconsistency in
	2018.		meta-analyses. Bmi. 2003 327(7414):557-60
22.	Finnish Arthroplasty Register. Total hip		DOI: 10.1136/bmj.327.7414.557.
	and knee arthroplasty report 2017.	33.	Viechtbauer W. Conducting meta-analyses
	Available from:	,,,	in R with metafor package. Journal of
	https://www2.thl.fi/endo/report/#data/kn		Statistical Software. 2010 36(3):1-48.
	ee topmodels. Accessed on July 14, 2018.	34.	Danish Knee Arthroplasty Register.
23.	The New Zealand Joint Registry. Eighteen	2.	Annual report 2019. Available from:
,	year report January 1999 to December		https://www.sundhed.dk/content/cms/99
	2016. Available from:		/4699_dkr-aarsrapport-2019_til-
	https://nzoa.org.nz/system/files/DH7827_		offentliggoerelse.pdf. Accessed on 9
	NZJR_2017_Report_v4_26Oct17.pdf.		December, 2019.
	Accessed on July 10, 2018.	35.	Dutch Arthroplasty Register (LROI).
24.	Norwegian Arthroplasty Register.		Online LROI annual report 2019 - PDF.
	Norwegian National Advisory Unit on		Available from: http://www.lroi-
	Arthroplasty and Hip Fractures: 2018		rapportage.nl/media/pdf/PDF%20Online
	Annual Report. Available from:		%20LROI%20annual%20report%202019.p
	http://nrlweb.ihelse.net/eng/Rapporter/R		df. Accessed on December 9, 2019.
	eport2018 english.pdf. Accessed on July	36.	Melvin T, Torre M. New medical device
	17, 2018.		regulations: the regulator's view. EFORT
25.	Swedish Knee Arthroplasty Register.		Open Rev. 2019 4(6):351-356 DOI:
	Annual Report 2017. Available from:		10.1302/2058-5241.4.180061.
	http://www.myknee.se/pdf/SVK_2017_Eng	37.	Faro LM, Huiskes R. Quality assurance of
	<u>1.0.pdf</u> . Accessed on July 14, 2018.		joint replacement. Legal regulation and
26.	National Joint Registry for England Wales		medical judgement. Acta Orthop Scand
	Northern Ireland and the Isle of Man. 14th		Suppl. 1992 250:1-33.
	Annual Report 2017. Available from:	38.	Hirst A, Philippou Y, Blazeby J, Campbell
	https://reports.njrcentre.org.uk/Portals/o/		B, Campbell M, Feinberg J, et al. No
	PDFdownloads/NJR%2014th%20Annual%		Surgical Innovation Without Evaluation:
	20Report%202017.pdf. Accessed on July		Evolution and Further Development of the
	20, 2018.		IDEAL Framework and Recommendations.
27.	Regional Register of Orthopaedic		Ann Surg. 2019 269(2):211-220 DOI:
-	Prosthetic Implantology. Hip, knee and		10.1097/sla.00000000002794.
	shoulder arthroplasty in Emilia-Romagna	39.	OECD (2017). Healthcare at a Glance 2017:
	Region (Italy) 2000-2016. Available from:		OECD Indicators. Available from:
			https://www.oecd-

ilibrary.org/docserver/health_glance-2017en.pdf?expires=1576672859&id=id&accna me=oido48188&checksum=4CFC37490313 AEE20A3B7DCAE71C4DC4. Accessed on August 23, 2019.

- Nieuwenhuijse MJ, Nelissen RG. [The introduction of new medical devices] [In Dutch]. Ned Tijdschr Geneeskd. 2015 159:A8652.
- Beyond Compliance Advisory Group. Beyond Compliance. Available from: www.beyondcompliance.org.uk/. Accessed on August 23, 2019.
- 42. Otten V, Maguire GQ, Jr., Noz ME, Zeleznik MP, Nilsson KG, Olivecrona H.

Are CT Scans a Satisfactory Substitute for the Follow-Up of RSA Migration Studies of Uncemented Cups? A Comparison of RSA Double Examinations and CT Datasets of 46 Total Hip Arthroplasties. Biomed Res Int. 2017 2017;3681458 DOI: 10.1155/2017/3681458.

Bozic KJ, Kurtz SM, Lau E, Ong K, Chiu V, Vail TP, et al. The epidemiology of revision total knee arthroplasty in the United States. Clin Orthop Relat Res. 2010 468(1):45-51 DOI: 10.1007/s11999-009-0945-0.

43.

Supplemental data

Non-RSA-Tested Designs at 5-year Follow-Up Revision rate with 95% confidence interval Australia ¹ Constantial Antonia (Alleges, Sydney, Anatonia) ² Constantial Antonials (Microsoft Orthopolitis, Arthopoliti, USA Sonstantial Antonials (Microsoft, Stechardine) Constantial Antonials (Microsoft Antonia), Tailingun, Garmany) Bodh⁴ Consent Assiste (Molece Internation), Canit I OME Prinary⁴ Consent American (Molece Internation) Gamin (Personnel CR: Molece American (Smith & Nepher) Mole⁴ Consent CR: Molec American (Smith & Nepher) Genein II Dasinet¹⁶ Concented DR Aneralis (Stelds & Nighew) Genein & Dasinet¹⁶ Concented Polyherinal (Stelds & Nighew) Josep Concent¹⁶ Concented National (Stelds & Nighew) Kannese Flat¹⁶ Cancented Aneralis (Steljster, Malwad, NJ, USA) USA¹⁶ Concented Concented (National Origins Malwad, NJ, USA) USA¹⁶ Concented Concented (National Origins (Stellster)) Commend Of Methic Assessing (Darky Speaker) Commend PE Methic Assessing (Darky Speaker) in IP*Commend OR Assessing (Smith & Nephow) in IP*Commend OR Assessing (Smith & Nephow) in IP*Commend PE Assessing (Smith & Nephow) ⁴⁴ Constituti Antibulis (Clement Hamet, Warsson, DK USA) of Kasa II¹⁴ Canonizal Antibulis (Zimmer Hamet) of Kasa II¹⁴ Canonizal Antibulis (Zimmer Hamet) of Kasa II¹⁴ Canonizal (El Antobelo (Clement Hamet)) Spatial-BBP Cancerd FS Martin, (Encode, Gamerik, BFCEpper Cancerd FS Math, Annula, (Daper Synthe) BFCEpper Cancerd FS Math, (Daper Synthe) Distinguist Cancerd FS Annual, (Daper Indust) Beth Concerd Annula, (Edds Debuges Indust) Beth Concerd Annula, (Edds Debuges In Torberg, SS Senser Concerd Annula, (Edds), (Mater, Freed Senser, WWT Cancerd C: Bachad, (Optica) Senser WWT Cancerd C: Bachad, (Optica) Sungar^a Consent PS Autoda (Stryker) Songar^a Consent PS Autoda (Stryker) Sens 1909^a Consent PS Autoda (Stryker) Nasir^a Consent (PS Autoda (Sensor Benct) Masir^a Consent PS Autoda (Sensor Benct) Masir^a Consent (PS Autoda (Sensor Benct) Autoralian¹⁰ Incommentar American (Magnet Advance¹⁰ Decemented American (Mergen Deleopelice) Advance¹⁰ Decemented American (Wright Medical Group Inc Calendra¹⁰ Uncommentar American (B Datas American) General (I¹⁰ Decementaria) (St. Barrelah (Decementaria) General (I¹⁰ Decementaria) (St. Barrelah (Decementaria)) (Decementaria) (Decementaria)) (Decementaria) (Decementaria) (Decementaria) (Decementaria) (Decementaria)) (Decementaria) (D Press, Annual A. Mark, A. Markan, Linner A. Pophery, General HT Unconstell (FS: Annualis (Statistic & Pophery) ECP Unconstell (FR: Annualis (DePpy Synthes) ECP Unconstell (FR: Molds: Annualis (DePpy Synthes) ECR Dec14* Unconstell (FR: Molds: Annualis (DePpy Synthes) ECR Dec14* Unconstell (FR: Molds: Annualis (DePpy Synthes) LCS Decla* Generated Melik Australia (Delvy Syste Macia* Decement Averatia (General Benet) Namel Rev II* Decement Averatio (General Benet) Namel Rev II* Decement Averatio (General Benet) Namel Rev II* Decement OK Averatio (General Benet) Name Rev* Decement OK Averatio (General Benet) Negar Per⁴ Decensiol CL Justich (Denner Bone) Negar Per VM CL Averale (Denner Bone) Negar UH Pho TMP PB Annola (Dinner Bone) 380 Dech⁴ Decensiol CR Annola (Deput) 380 Dech⁴ Decensiol CR Moldi, Annola (Deput Soulo) 297 Span MEP⁴ Decensiol CR Moldi, Annola (Deput Soulo) 297 Span MEP⁴ Decensiol CR Moldi, Annola (Deput Soulo) PECSpea MIT Deals* December (Childe Aurola (Childe Inner (909° Uncommed CE Antolas (Brylor) Saria (909° Uncommed Pi-Intrafa (Brylor) Seria (909° Uncommed Pi-Intrafa (Brylor) Tsahlar⁴ Uncommed PS Antolas (Strylor) Matte⁴ Uncommed CE Antolas (Strylor) Matte⁴ Uncommed CE Antolas (Strylor) Sugarrer, Commune C. Sarrela (Comm. Rent) Vargaret, Commune C. Sarrela (Comm. Rent) Activities, "Hybrid Astrola (Mary): Advance IP Blood Astrola (Mary): Gammin IP Blood (C. Astrola (Conth & Supher)) ILCP Hold DK Austalia (Editory Syndau ILCP Hold DKA Austalia (Editory Syndau) ILCE Inter[®] Hilterd Methic Assestis (Editory Syndau) ILCE Inter[®] Hilterd Methic Assestis (Editory Syndau) Concess (FF Hold CK Austalia (Editory Syndau) Material Methic Assestis (Editory Benefic) Material Methic Assestis (Editory Benefic) Material Methic Assestis (Editory Benefic) Material Methica (Editory Benefic) Material Methica (Editory Benefic) Material Methica (Editory Benefic) Material Methica (Editory Benefic) Mich Synda[®] Mitt (Hold CK Assestis (Editory Benefic) Mich Synda[®] Mitt (Hold CK Assestis (Editory Benefic) Mich Synda[®] Mitt (Hold CK Assestis (Editory Benefic) PEC Signal⁶ MEE Daulis Holmil PS Annulas (Depur Sprit Paulis⁶ Hybrid Holmin Annulas (Daulis & Naphow) ERIC⁶ Hybrid Annulas (Education Companie Tachandegy) Samia (1994) Hybrid (D. Annulas Olitykov) Samia (2009) Hybrid (D. Annulas Olitykov) Samia (2009) Hybrid (D. Annulas Olitykov) Baran 2009" Hotel CR Anenda (Dryker) Sergint[®] Hotel PS Anenda (Dryker) Saria 2009" Hotel PS Anenda (Sryker) Made" Hotel CR Anenda (Drener Brend) Made" Hotel CR Anenda (Drener Brend) AGC Dual Arisolar Common Hendel AGC Dual Arisolar Common Federal (Jamme Home) AGC '97 Consued Fuland (Jamme Home) AGC '97 Moniel" Consolid Falad (Jamme Home) Assess?" Consult Falad (Style) Assess?" Consult Falad (Style) Finland Age⁴⁴ Concept Field Conference Only, Rockeds, Com Genesis II* Concept (C Fisland (Sould: & Nephere) Genesis II* Concept (C Fisland (Sould: & Nephere) Kasenate: Fall* Concept Fisland (Soviet) Match* Concept Fisland (Soviet) Match* Concept Fisland (Soviet) Source Fis: TWP (C) Fisland (Soviet) Source Fis: TWP (C) Fisland (Soviet) Negar Fox 199° CE Frand (Denie) Bond) Nogar DM* PSFiled (Denie) Roma) RCI Inters DM* Consold Field (Dolle) RCI Inters DM* Consold Field (Dolle) RCI Span Stabilish* Consold PS Field (Dolle System) RCI Span Stabilish* Consold PS Field (Doll System) RCI Span Stabilish* Consold PS Field (Doll System) TON Concession Fields" TON Concession Fields (Concession) Tradition" Concession 715 Fields (Concession) Varguards" (Traditional PS Fields) (Concessions) Formation Forwards (Fields) Fields) General II High Honort⁴ FS Incite Romagne (Smith & Ne Varguad⁴ FS Enda Romagne (Sannor Hender) General MC 3⁴⁴ Tendin Romagne (Link, Henders, Germany TC File SD Enlanz⁴ Tendin Romagne (Smith & Higher) Emilia-Romagna (Italy) Parlis Candorsing¹⁰ Darlis Romagna (South-& Naphor) PCC R¹⁰ 15 Mobile Davids Tamagna (Dopt Syndoc) Gaussis II¹⁰ CR Sindia Parsagna (Doub & Naghors) Gauss II¹⁰ CR Sindia Romagna (Malis Obto, Comman, Inde) Vangand¹⁰ Layerd CR Davids Farsagna (Dower Bornet) Hand Kinda Komman¹⁰ Roughls" Daille Romanne (Dirit Group, Classeveter, UR) PTC*P5 Edit Romanne (Darry Syndro) Group TeXXT* Roman (Darry Syndro) Oracle TeXXX*P Romanne (Aster Michael, Miles, Indy Profil**95 Daille Romanne (Aster Michael, Miles, Indy Advance** Nobile/Proct Darlies Romagne (Microsoft Ochowski Advance** Nobile/Proct Darlies Romagne (Microsoft Ochowski Strenger) . Multi-resolution Holing and Khangar Ulthaneo Changa Maran Andrea San Andrea San Angel Ulthaneo Changa Maran Angel Angel Angel Angel Angel Angel Maran Maran Angel Angel Angel Maran Maran Angel Angel Angel Angel Angel Angel Angel Maran Maran Angel Angel Maran Maran Maran Angel Maran Maran Maran Maran Angel Maran Nongen Unconnected New Zealand, Communities LCS Community New Zealand, Delving Spotkop LCS Holed New Zealand (Delving Spotkop LCS Holed New Zealand (Delving System) ECS Homowards New Zealand (Delving Spotkop) Scorege Community New Zealand (Delving Spotkop) New Zealand + -----EX Constant New Zoland Learner News) Scorps Constant New Zoland (Zerrer News) Isodification Constant/New Zoland (Zerrer News) Adusc# Constant New Zoland (Mersyan Orbopolici) Adusc# Constant New Zoland (Mersyan Orbopolici) FrancesCan MEP Synder (Inder Orthque FC Signe* APT Sender (Dans Synder) FC NeM* Sender (Dans Sender) Sweden PEC Mobile[®] Sweden (Dasse Studiet) PEC Sprat[®] Connected OD Hobbis UK (Dasse Synthes) Makawa M[®] Connected OF Hoad UK (Merupent Orthopades) Makawa PH[®] Connected OF Hoad UK (Diffuse Ancodag) Similari[®] Connected OF Hoad UK (Diffuse Ancodag) Sweden Showeljan Kon[®] Connected OF Hobbis UK (Diffuse Ancodag) United Kingdom Hermighter Kanet¹⁰ Connected CE Michile UK (2) Brauer, Bacondydar Konet¹⁰ Holdel CD Michile UK (2) Brauer, A 19 Connected (20 Front UK (Smith & Hopbert) 10 Oannet¹⁰ Connected (20 Front UK (Smith & Hopbert) 10 Oannet¹⁰ Connected (20 Front UK (Smith & Hopbert) 10 Oannet¹⁰ Connected (20 Front M (Smith & Hopbert) Kinemas* Commod CH Final UK (Stryles) LCP* Hybrid CH Mick UK (Brilly: Soukes) LCS Complex* Consoled CH Mickie UK (Brilly Symbol) LCS Complex* Consoled CH Mick UK (Drilly Symbol) MR* Consoled CH Fixed UK (Tables) Collegenders) Marine¹⁰ Consented CM Fixed UK (Zeniner Barrier) MR2⁴⁰ Consented PS Fixed (Seniner) Natural Kare IP Consented CR Fixed UK (Seniner Bar PC Signa⁴⁰ Consented PS Matche UK (Dapas Tyrelos NY Tirgun⁴⁰ Consented PS Matche UK (Dapas Tyrelos Scorpie⁴ Consent CH Mildle UK (Strylor) Scorpie⁴ Consent PS Fault UK (Strylor) Scorpie⁴ Consent PS Malit: UK(Strylor) Scorpie⁴ Bylen CH Faul UK (Strylor) TC Pie⁴ Consent CK Faul UK (Strylor) TC Pie⁴ Consent CK Faul UK (Strylor) -TC Plus⁴ Centrated CR Multie UK (Double & Nepher TC Plus⁴ Isyland CR Multie UK (Enable & Nephere) Vanguard⁴ Constant PS Fixed UK (Zentrat Borred Schurze MP States⁴ Constant CR (Zentrat Borred) Schurze MP States⁴ Constant CR (Plust UK (Mersy ... 3.6% $I^2 = 96\%$ Pooled revision rate of non-RSA tested designs 0 2 4 6 8 10 12 14 Revision rate at 5-year follow-up in percentages

Supplemental II.I

This forest plot shows revision rates with 95% CIs at 5 years of follow-up of the non-RSA-tested TKA designs, subdivided per registry, and the pooled revision rate at 5 years of follow-up with 95% CI. The heterogeneity between TKA designs stands out, with an I² of 96%.



Supplemental II.II

Forest plot showing the revision rates with 95% confidence intervals at 5-year follow-up of the RSAtested TKP designs subdivided per registry, and the pooled revision rate at 5-year follow-up with 95% confidence interval. The heterogeneity between the TKP designs stands out with a 1² of 98%.

Non-RSA-Tested Designs at 10-year Follow-Up



Supplemental II.III

Forest plot showing the revision rates with 95% confidence intervals at 10-year follow-up of the non-RSA-tested TKP designs subdivided per registry, and the pooled revision rate at 10-year follow-up with 95% confidence interval. The heterogeneity between TKP designs stands out with a l² of 97%



Supplemental II.IV

Forest plot showing the revision rates with 95% confidence intervals at 10-year follow-up of the RSAtested TKP designs subdivided per registry, and the pooled revision rate at 10-year follow-up with 95% confidence interval. The heterogeneity between the TKP designs stands out with a 1² of 97%.